

## Product datasheet for **TL308863**

### TFIP11 Human shRNA Plasmid Kit (Locus ID 24144)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | TFIP11 Human shRNA Plasmid Kit (Locus ID 24144)  |
| Locus ID:                 | 24144  |
| Synonyms:                 | bK445C9.6; NTR1; Spp382; STIP; STIP-1; TIP39   |
| Vector:                   | pGFP-C-shLenti (TR30023)   |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)   |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Lentiviral plasmids  |
| Components:               | TFIP11 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 24144).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.  |
| RefSeq:                   | <a href="#">NM_001008697</a> , <a href="#">NM_012143</a> , <a href="#">NM_001346857</a> , <a href="#">NM_001346858</a> , <a href="#">NM_001346859</a> ,<br><a href="#">NM_001346861</a> , <a href="#">NM_001346862</a> , <a href="#">NM_001008697.1</a> , <a href="#">NM_001008697.2</a> , <a href="#">NM_012143.1</a> ,<br><a href="#">NM_012143.2</a> , <a href="#">NM_012143.3</a> , <a href="#">BC011599</a> , <a href="#">BC033080</a> , <a href="#">BC033080.1</a> , <a href="#">NM_012143.4</a>   |
| UniProt ID:               | <a href="#">Q9UBB9</a>   |
| Summary:                  | This gene encodes a protein component of the spliceosome that promotes the release of the lariat-intron during late-stage splicing through the recruitment of a pre-mRNA splicing factor called DEAH-box helicase 15. The encoded protein contains a G-patch domain, a hallmark of RNA-processing proteins, that binds DEAH-box helicase 15. This protein contains an atypical nuclear localization sequence as well as a nuclear speckle-targeting sequence, enabling it to localize to distinct speckled regions within the cell nucleus. Polymorphisms in this gene are associated with dental caries suggesting a role in amelogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2016] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).