

## Product datasheet for **TL308865**

### Spasmolytic Polypeptide (TFF2) Human shRNA Plasmid Kit (Locus ID 7032)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Spasmolytic Polypeptide (TFF2) Human shRNA Plasmid Kit (Locus ID 7032)
Locus ID:	7032
Synonyms:	SML1; SP
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TFF2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 7032). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_005423</a> , <a href="#">NM_005423.1</a> , <a href="#">NM_005423.2</a> , <a href="#">NM_005423.3</a> , <a href="#">NM_005423.4</a> , <a href="#">BC032820</a> , <a href="#">BC032820.1</a> , <a href="#">BM817851</a> , <a href="#">BM853946</a> , <a href="#">NM_005423.5</a>
UniProt ID:	<a href="#">Q03403</a>
Summary:	Members of the trefoil family are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain that contains three conserved disulfides. They are stable secretory proteins expressed in gastrointestinal mucosa. Their functions are not defined, but they may protect the mucosa from insults, stabilize the mucus layer and affect healing of the epithelium. The encoded protein inhibits gastric acid secretion. This gene and two other related trefoil family member genes are found in a cluster on chromosome 21. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).