

Product datasheet for **TL308893V**

Cripto1 (TDGF1) Human shRNA Lentiviral Particle (Locus ID 6997)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cripto1 (TDGF1) Human shRNA Lentiviral Particle (Locus ID 6997)
Locus ID:	6997
Synonyms:	CR; CR-1; CRGF; CRIPTO
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	TDGF1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC022393 , NM_001174136 , NM_003212 , NM_003212.1 , NM_003212.2 , NM_003212.3 , NM_001174136.1 , BC067844 , BC067844.1 , NM_001174136.2 , NM_003212.4
UniProt ID:	P13385
Summary:	This gene encodes an epidermal growth factor-related protein that contains a cripto, FRL-1, and cryptic domain. The encoded protein is an extracellular, membrane-bound signaling protein that plays an essential role in embryonic development and tumor growth. Mutations in this gene are associated with forebrain defects. Pseudogenes of this gene are found on chromosomes 2, 3, 6, 8, 19 and X. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).