

## Product datasheet for **TL308934**

### **TBC1D5 Human shRNA Plasmid Kit (Locus ID 9779)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	TBC1D5 Human shRNA Plasmid Kit (Locus ID 9779)
Locus ID:	9779
Synonyms:	KIAA0210
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TBC1D5 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9779). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001134380</a> , <a href="#">NM_001134381</a> , <a href="#">NM_014744</a> , <a href="#">NM_001349073</a> , <a href="#">NM_001349074</a> , <a href="#">NM_001349075</a> , <a href="#">NM_001349076</a> , <a href="#">NM_001349077</a> , <a href="#">NM_001349078</a> , <a href="#">NM_001349079</a> , <a href="#">NM_001349080</a> , <a href="#">NM_001349081</a> , <a href="#">NM_001349082</a> , <a href="#">NM_001349083</a> , <a href="#">NM_001349084</a> , <a href="#">NM_001349085</a> , <a href="#">NM_001349086</a> , <a href="#">NM_001349087</a> , <a href="#">NM_001349088</a> , <a href="#">NM_001349089</a> , <a href="#">NM_001349090</a> , <a href="#">NM_001349091</a> , <a href="#">NM_014744.1</a> , <a href="#">NM_014744.2</a> , <a href="#">NM_001134380.1</a> , <a href="#">NM_001134381.1</a> , <a href="#">BC013145</a> , <a href="#">BM976960</a> , <a href="#">NM_001134381.2</a>
UniProt ID:	<a href="#">Q92609</a>
Summary:	May act as a GTPase-activating protein (GAP) for Rab family protein(s). May act as a GAP for RAB7A. Can displace RAB7A and retromer CSC subcomplex from the endosomal membrane to the cytosol; at least retromer displacement seems to require its catalytic activity (PubMed:19531583, PubMed:20923837). Required for retrograde transport of cargo proteins from endosomes to the trans-Golgi network (TGN); the function seems to require its catalytic activity. Involved in regulation of autophagy (PubMed:22354992). May act as a molecular switch between endosomal and autophagosomal transport and is involved in reprogramming vesicle trafficking upon autophagy induction. Involved in the trafficking of ATG9A upon activation of autophagy. May regulate the recruitment of ATG9A-AP2-containing vesicles to autophagic membranes (PubMed:24603492).[UniProtKB/Swiss-Prot Function]


[View online »](#)

<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .
<b>Performance Guaranteed:</b>	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>