

Product datasheet for **TL308987**

hnRNP Q (SYNCRIP) Human shRNA Plasmid Kit (Locus ID 10492)

Product data:

Product Type:	shRNA Plasmids
Product Name:	hnRNP Q (SYNCRIP) Human shRNA Plasmid Kit (Locus ID 10492)
Locus ID:	10492
Synonyms:	GRY-RBP; GRYRBP; hnRNP-Q; HNRNPQ; HNRNPQ1; NSAP1; PP68
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SYNCRIP - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10492). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001159673 , NM_001159674 , NM_001159675 , NM_001159676 , NM_001159677 , NM_001253771 , NM_006372 , NM_006372.1 , NM_006372.2 , NM_006372.3 , NM_006372.4 , NM_001159673.1 , NM_001159676.1 , NM_001159677.1 , NM_001159675.1 , NM_001159674.1 , NM_001253771.1 , BC032643 , BC009176 , BC015575 , BC021932 , BC024283 , BC040844 , NM_001159674.2 , NM_001253771.2 , NM_001159677.2 , NM_006372.5
UniProt ID:	O60506
Summary:	This gene encodes a member of the cellular heterogeneous nuclear ribonucleoprotein (hnRNP) family. hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA (hnRNA) and regulate alternative splicing, polyadenylation, and other aspects of mRNA metabolism and transport. The encoded protein plays a role in multiple aspects of mRNA maturation and is associated with several multiprotein complexes including the apoB RNA editing-complex and survival of motor neurons (SMN) complex. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the short arm of chromosome 20. [provided by RefSeq, Dec 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).