

Product datasheet for **TL309044**

STOM Human shRNA Plasmid Kit (Locus ID 2040)

Product data:

Product Type:	shRNA Plasmids
Product Name:	STOM Human shRNA Plasmid Kit (Locus ID 2040)
Locus ID:	2040
Synonyms:	BND7; EPB7; EPB72
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	STOM - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 2040). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001270526 , NM_001270527 , NM_004099 , NM_198194 , NR_073037 , NM_004099.1 , NM_004099.2 , NM_004099.3 , NM_004099.4 , NM_004099.5 , NM_198194.1 , NM_198194.2 , NM_001270527.1 , NM_001270526.1 , BC010703 , BC010703.1 , BM451470 , BM925356 , NM_001270526.2 , NM_198194.3
UniProt ID:	P27105
Summary:	This gene encodes a member of a highly conserved family of integral membrane proteins. The encoded protein localizes to the cell membrane of red blood cells and other cell types, where it may regulate ion channels and transporters. Loss of localization of the encoded protein is associated with hereditary stomatocytosis, a form of hemolytic anemia. There is a pseudogene for this gene on chromosome 6. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).