

## Product datasheet for **TL309117**

### **SREBP1 (SREBF1) Human shRNA Plasmid Kit (Locus ID 6720)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	SREBP1 (SREBF1) Human shRNA Plasmid Kit (Locus ID 6720)
<b>Locus ID:</b>	6720
<b>Synonyms:</b>	bHLHd1; HMD; IFAP2; SREBP1
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	SREBF1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6720). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<a href="#">NM_001005291</a> , <a href="#">NM_001321096</a> , <a href="#">NM_004176</a> , <a href="#">NM_004176.1</a> , <a href="#">NM_004176.2</a> , <a href="#">NM_004176.3</a> , <a href="#">NM_004176.4</a> , <a href="#">NM_001005291.1</a> , <a href="#">NM_001005291.2</a> , <a href="#">BC057388</a> , <a href="#">BC057388.1</a> , <a href="#">BC063281</a> , <a href="#">BC063281.1</a> , <a href="#">BC023621</a> , <a href="#">BC026962</a> , <a href="#">NM_001005291.3</a>
<b>UniProt ID:</b>	<a href="#">P36956</a>
<b>Summary:</b>	This gene encodes a basic helix-loop-helix-leucine zipper (bHLH-Zip) transcription factor that binds to the sterol regulatory element-1 (SRE1), which is a motif that is found in the promoter of the low density lipoprotein receptor gene and other genes involved in sterol biosynthesis. The encoded protein is synthesized as a precursor that is initially attached to the nuclear membrane and endoplasmic reticulum. Following cleavage, the mature protein translocates to the nucleus and activates transcription. This cleavage is inhibited by sterols. This gene is located within the Smith-Magenis syndrome region on chromosome 17. Alternative promoter usage and splicing result in multiple transcript variants, including SREBP-1a and SREBP-1c, which correspond to RefSeq transcript variants 2 and 3, respectively. [provided by RefSeq, Nov 2017]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).