

Product datasheet for **TL309184**

SORCS1 Human shRNA Plasmid Kit (Locus ID 114815)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SORCS1 Human shRNA Plasmid Kit (Locus ID 114815)
Locus ID:	114815
Synonyms:	hSorCS
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SORCS1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 114815). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001013031 , NM_001206569 , NM_001206570 , NM_001206571 , NM_001206572 , NM_052918 , NM_052918.1 , NM_052918.2 , NM_052918.3 , NM_052918.4 , NM_001013031.1 , NM_001013031.2 , NM_001206570.1 , NM_001206571.1 , NM_001206569.1 , NM_001206572.1 , BC131597 , NM_001013031.3 , NM_001206569.2 , NM_052918.5 , NM_001206570.2 , NM_001206571.2 , NM_001206572.2
UniProt ID:	Q8WY21
Summary:	This gene encodes one family member of vacuolar protein sorting 10 (VPS10) domain-containing receptor proteins. The VPS10 domain name comes from the yeast carboxypeptidase Y sorting receptor Vps10 protein. Members of this gene family are large with many exons but the CDS lengths are usually less than 3700 nt. Very large introns typically separate the exons encoding the VPS10 domain; the remaining exons are separated by much smaller-sized introns. These genes are strongly expressed in the central nervous system. Two of the five family members (sortilin and sortilin-related receptor) are synthesized as preproteins; it is not yet known if this encoded protein is also a preproprotein. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]



[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>