

Product datasheet for **TL309239**

SMC4 Human shRNA Plasmid Kit (Locus ID 10051)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SMC4 Human shRNA Plasmid Kit (Locus ID 10051)
Locus ID:	10051
Synonyms:	CAPC; CAPC, SMC4L1, hCAP-C; chromosome-associated polypeptide C; hCAP-C; SMC4 structural maintenance of chromosomes 4-like 1; SMC4L1; structural maintenance of chromosomes 4
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SMC4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10051). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001002799 , NM_001002800 , NM_001288753 , NM_005496 , NM_001002800.1 , NM_001002800.2 , NM_005496.1 , NM_005496.2 , NM_005496.3 , NM_001288753.1 , NM_001002799.1 , BC001557 , BC065259 , BC070161 , BC106033 , BC146370 , BC148796 , NM_001288753.2 , NM_001002800.3
UniProt ID:	Q9NTJ3
Summary:	This gene belongs to the 'structural maintenance of chromosomes' (SMC) gene family. Members of this gene family play a role in two changes in chromosome structure during mitotic segregation of chromosomes- chromosome condensation and sister chromatid cohesion. The protein encoded by this gene is likely a subunit of the 13S condensin complex, which is involved in chromosome condensation. A pseudogene related to this gene is located on chromosome 2. [provided by RefSeq, Jun 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).