

## Product datasheet for **TL309295**

### SLC6A14 Human shRNA Plasmid Kit (Locus ID 11254)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SLC6A14 Human shRNA Plasmid Kit (Locus ID 11254)
Locus ID:	11254
Synonyms:	BMIQ11
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SLC6A14 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 11254). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_007231</a> , <a href="#">NM_007231.1</a> , <a href="#">NM_007231.2</a> , <a href="#">NM_007231.3</a> , <a href="#">NM_007231.4</a> , <a href="#">BC093710</a> , <a href="#">BC093710.1</a> , <a href="#">BC093712</a> , <a href="#">NM_007231.5</a>
UniProt ID:	<a href="#">Q9UN76</a>
Summary:	This gene encodes a member of the solute carrier family 6. Members of this family are sodium and chloride dependent neurotransmitter transporters. The encoded protein transports both neutral and cationic amino acids. This protein may also function as a beta-alanine carrier. Mutations in this gene may be associated with X-linked obesity. A pseudogene of this gene is found on chromosome X.[provided by RefSeq, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL309295A virus into HEK293 cells. TL309295A virus was prepared using lenti-shRNA TL309295A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL309295B virus into HEK293 cells. TL309295B virus was prepared using lenti-shRNA TL309295B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL309295C] virus into HEK293 cells. [TL309295C] virus was prepared using lenti-shRNA [TL309295C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL309295D] virus into HEK293 cells. [TL309295D] virus was prepared using lenti-shRNA [TL309295D] and [TR30037] packaging kit.