

Product datasheet for **TL309363**

SLC25A25 Human shRNA Plasmid Kit (Locus ID 114789)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SLC25A25 Human shRNA Plasmid Kit (Locus ID 114789)
Locus ID:	114789
Synonyms:	MCSC; PCSCL; SCAMC-2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SLC25A25 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 114789). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001006641 , NM_001006642 , NM_001006643 , NM_001265614 , NM_001330988 , NM_052901 , NR_049766 , NM_052901.1 , NM_052901.2 , NM_052901.3 , NM_052901.4 , NM_001006641.1 , NM_001006641.2 , NM_001006641.3 , NM_001006642.1 , NM_001006642.2 , NM_001006642.3 , NM_001265614.1 , NM_001265614.2 , NM_001006643.1 , BC089448 , BC005163 , BC036889 , BC044631 , BC061925 , BC103930 , BC103931 , BC103932 , BC103933 , BM689843
UniProt ID:	Q6KCM7
Summary:	The protein encoded by this gene belongs to the family of calcium-binding mitochondrial carriers, with a characteristic mitochondrial carrier domain at the C-terminus. These proteins are found in the inner membranes of mitochondria, and function as transport proteins. They shuttle metabolites, nucleotides and cofactors through the mitochondrial membrane and thereby connect and/or regulate cytoplasm and matrix functions. This protein may function as an ATP-Mg/Pi carrier that mediates the transport of Mg-ATP in exchange for phosphate, and likely responsible for the net uptake or efflux of adenine nucleotides into or from the mitochondria. Alternatively spliced transcript variants encoding different isoforms with a common C-terminus but variable N-termini have been described for this gene. [provided by RefSeq, Jul 2012]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).