

## Product datasheet for **TL309372V**

### SLC22A8 Human shRNA Lentiviral Particle (Locus ID 9376)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	SLC22A8 Human shRNA Lentiviral Particle (Locus ID 9376)
Locus ID:	9376
Synonyms:	OAT3
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	SLC22A8 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001184732</a> , <a href="#">NM_001184733</a> , <a href="#">NM_001184736</a> , <a href="#">NM_004254</a> , <a href="#">NM_004254.1</a> , <a href="#">NM_004254.2</a> , <a href="#">NM_004254.3</a> , <a href="#">NM_001184736.1</a> , <a href="#">NM_001184733.1</a> , <a href="#">NM_001184732.1</a> , <a href="#">BC022387</a> , <a href="#">BC022387.1</a> , <a href="#">NM_004254.4</a>
UniProt ID:	<a href="#">Q8TCC7</a>
Summary:	This gene encodes a protein involved in the sodium-independent transport and excretion of organic anions, some of which are potentially toxic. The encoded protein is an integral membrane protein and appears to be localized to the basolateral membrane of the kidney. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene. [provided by RefSeq, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).