

## Product datasheet for **TL309436**

### SIPA1L3 Human shRNA Plasmid Kit (Locus ID 23094)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SIPA1L3 Human shRNA Plasmid Kit (Locus ID 23094)
Locus ID:	23094
Synonyms:	CTRCT45; SPAL3; SPAR3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SIPA1L3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 23094). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_015073</a> , <a href="#">NM_015073.1</a> , <a href="#">NM_015073.2</a> , <a href="#">BC036863</a> , <a href="#">BC150620</a> , <a href="#">BM193478</a>
UniProt ID:	<a href="#">O60292</a>
Summary:	This gene belongs to the signal induced proliferation associated 1 family of genes, which encode GTPase-activating proteins specific for the GTP-binding protein Rap1. Rap1 has been implicated in regulation of cell adhesion, cell polarity, and organization of the cytoskeleton. Like other members of the family, the protein encoded by this gene contains RapGAP and PDZ domains. In addition, this protein contains a C-terminal leucine zipper domain. This gene is proposed to function in epithelial cell morphogenesis and establishment or maintenance of polarity. Consistently, expression of the protein in cell culture showed localization to cell-cell borders in apical regions, and downregulation of the gene in 3D Caco2 cell culture resulted in abnormal cell polarity and morphogenesis. Allelic variants of this gene have been associated with congenital cataracts in humans. [provided by RefSeq, Feb 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).