

## Product datasheet for **TL309480**

### SFTPC Human shRNA Plasmid Kit (Locus ID 6440)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SFTPC Human shRNA Plasmid Kit (Locus ID 6440)
Locus ID:	6440
Synonyms:	BRICD6; PSP-C; SFTP2; SMDP2; SP-C; SP5
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SFTPC - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6440). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001172357</a> , <a href="#">NM_001172410</a> , <a href="#">NM_001317778</a> , <a href="#">NM_001317779</a> , <a href="#">NM_001317780</a> , <a href="#">NM_003018</a> , <a href="#">NM_003018.1</a> , <a href="#">NM_003018.2</a> , <a href="#">NM_003018.3</a> , <a href="#">NM_001172357.1</a> , <a href="#">NM_001172410.1</a> , <a href="#">BC005913</a> , <a href="#">BC005913.1</a> , <a href="#">BC105594</a> , <a href="#">BM993053</a> , <a href="#">NM_001172357.2</a> , <a href="#">NM_001172410.2</a>
UniProt ID:	<a href="#">P11686</a>
Summary:	This gene encodes the pulmonary-associated surfactant protein C (SPC), an extremely hydrophobic surfactant protein essential for lung function and homeostasis after birth. Pulmonary surfactant is a surface-active lipoprotein complex composed of 90% lipids and 10% proteins which include plasma proteins and apolipoproteins SPA, SPB, SPC and SPD. The surfactant is secreted by the alveolar cells of the lung and maintains the stability of pulmonary tissue by reducing the surface tension of fluids that coat the lung. Multiple mutations in this gene have been identified, which cause pulmonary surfactant metabolism dysfunction type 2, also called pulmonary alveolar proteinosis due to surfactant protein C deficiency, and are associated with interstitial lung disease in older infants, children, and adults. Alternatively spliced transcript variants encoding different protein isoforms have been identified.[provided by RefSeq, Feb 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).