

Product datasheet for **TL309489V**

SFRS3 (SRSF3) Human shRNA Lentiviral Particle (Locus ID 6428)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	SFRS3 (SRSF3) Human shRNA Lentiviral Particle (Locus ID 6428)
Locus ID:	6428
Synonyms:	SFRS3; SRp20
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	SRSF3 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_003017 , NR_036610 , NM_003017.1 , NM_003017.2 , NM_003017.3 , NM_003017.4 , BC069018 , BC069018.1 , BC000914 , BM833102 , NM_003017.5
UniProt ID:	P84103
Summary:	The protein encoded by this gene is a member of the serine/arginine (SR)-rich family of pre-mRNA splicing factors, which constitute part of the spliceosome. Each of these factors contains an RNA recognition motif (RRM) for binding RNA and an RS domain for binding other proteins. The RS domain is rich in serine and arginine residues and facilitates interaction between different SR splicing factors. In addition to being critical for mRNA splicing, the SR proteins have also been shown to be involved in mRNA export from the nucleus and in translation. Two transcript variants, one protein-coding and the other non-coding, have been found for this gene. [provided by RefSeq, Sep 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL309489A virus into HEK293 cells. TL309489A virus was prepared using lenti-shRNA TL309489A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL309489B virus into HEK293 cells. TL309489B virus was prepared using lenti-shRNA TL309489B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL309489C] virus into HEK293 cells. [TL309489C] virus was prepared using lenti-shRNA [TL309489C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL309489D] virus into HEK293 cells. [TL309489D] virus was prepared using lenti-shRNA [TL309489D] and [TR30037] packaging kit.