

## Product datasheet for **TL309491**

### SC35 (SRSF2) Human shRNA Plasmid Kit (Locus ID 6427)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SC35 (SRSF2) Human shRNA Plasmid Kit (Locus ID 6427)
Locus ID:	6427
Synonyms:	PR264; SC-35; SC35; SFRS2; SFRS2A; SRp30b
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SRSF2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6427). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001195427</a> , <a href="#">NM_003016</a> , <a href="#">NR_036608</a> , <a href="#">NM_003016.1</a> , <a href="#">NM_003016.2</a> , <a href="#">NM_003016.3</a> , <a href="#">NM_003016.4</a> , <a href="#">NM_001195427.1</a> , <a href="#">BC070086</a> , <a href="#">BC000339</a> , <a href="#">BC001303</a> , <a href="#">BC063570</a> , <a href="#">BC066958</a>
UniProt ID:	<a href="#">Q01130</a>
Summary:	The protein encoded by this gene is a member of the serine/arginine (SR)-rich family of pre-mRNA splicing factors, which constitute part of the spliceosome. Each of these factors contains an RNA recognition motif (RRM) for binding RNA and an RS domain for binding other proteins. The RS domain is rich in serine and arginine residues and facilitates interaction between different SR splicing factors. In addition to being critical for mRNA splicing, the SR proteins have also been shown to be involved in mRNA export from the nucleus and in translation. Two transcript variants encoding the same protein and one non-coding transcript variant have been found for this gene. In addition, a pseudogene of this gene has been found on chromosome 11. [provided by RefSeq, Sep 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).