

## Product datasheet for **TL309594**

### Syntenin (SDCBP) Human shRNA Plasmid Kit (Locus ID 6386)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Syntenin (SDCBP) Human shRNA Plasmid Kit (Locus ID 6386)
Locus ID:	6386
Synonyms:	MDA-9; MDA9; ST1; SYCL; TACIP18
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SDCBP - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6386). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001007067</a> , <a href="#">NM_001007068</a> , <a href="#">NM_001007069</a> , <a href="#">NM_001007070</a> , <a href="#">NM_005625</a> , <a href="#">NM_001330537</a> , <a href="#">NM_001348339</a> , <a href="#">NM_001348340</a> , <a href="#">NM_001348341</a> , <a href="#">NM_001007070.1</a> , <a href="#">NM_001007067.1</a> , <a href="#">NM_005625.1</a> , <a href="#">NM_005625.2</a> , <a href="#">NM_005625.3</a> , <a href="#">NM_001007068.1</a> , <a href="#">NM_001007069.1</a> , <a href="#">BC113674</a> , <a href="#">BC013254</a> , <a href="#">BC113676</a> , <a href="#">BC143915</a> , <a href="#">BC143916</a> , <a href="#">NM_005625.4</a>
UniProt ID:	<a href="#">O00560</a>
Summary:	The protein encoded by this gene was initially identified as a molecule linking syndecan-mediated signaling to the cytoskeleton. The syntenin protein contains tandemly repeated PDZ domains that bind the cytoplasmic, C-terminal domains of a variety of transmembrane proteins. This protein may also affect cytoskeletal-membrane organization, cell adhesion, protein trafficking, and the activation of transcription factors. The protein is primarily localized to membrane-associated adherens junctions and focal adhesions but is also found at the endoplasmic reticulum and nucleus. Alternative splicing results in multiple transcript variants encoding different isoforms. Related pseudogenes have been identified on multiple chromosomes. [provided by RefSeq, Jan 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).