

## Product datasheet for **TL309753**

### RPL17 Human shRNA Plasmid Kit (Locus ID 6139)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	RPL17 Human shRNA Plasmid Kit (Locus ID 6139)
Locus ID:	6139
Synonyms:	L17; PD-1; RPL23
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RPL17 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6139). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_000985</a> , <a href="#">NM_001035006</a> , <a href="#">NM_001199340</a> , <a href="#">NM_001199341</a> , <a href="#">NM_001199342</a> , <a href="#">NM_001199343</a> , <a href="#">NM_001199344</a> , <a href="#">NM_001199345</a> , <a href="#">NM_000985.1</a> , <a href="#">NM_000985.2</a> , <a href="#">NM_000985.3</a> , <a href="#">NM_000985.4</a> , <a href="#">NM_001035006.1</a> , <a href="#">NM_001035006.2</a> , <a href="#">NM_001199345.1</a> , <a href="#">NM_001199340.1</a> , <a href="#">NM_001199342.1</a> , <a href="#">NM_001199343.1</a> , <a href="#">NM_001199344.1</a> , <a href="#">BC000502</a> , <a href="#">BC000502.2</a> , <a href="#">BC017831</a> , <a href="#">BC017831.1</a> , <a href="#">BC066323</a> , <a href="#">BC066324</a> , <a href="#">BC106031</a> , <a href="#">NM_001369555</a> , <a href="#">NM_001369558</a> , <a href="#">NM_001369561</a> , <a href="#">NM_001369562</a> , <a href="#">NM_001369556</a> , <a href="#">NM_001369557</a> , <a href="#">NM_001369560</a> , <a href="#">NM_001369563</a> , <a href="#">NM_001199344.2</a> , <a href="#">NM_001199343.2</a> , <a href="#">NM_001199342.3</a> , <a href="#">NM_001199340.2</a> , <a href="#">NM_001199345.2</a> , <a href="#">NM_001199341.2</a> , <a href="#">NM_001035006.3</a> , <a href="#">NM_000985.5</a>
UniProt ID:	<a href="#">P18621</a>



[View online »](#)

<b>Summary:</b>	<p>Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L22P family of ribosomal proteins. It is located in the cytoplasm. This gene has been referred to as rpl23 because the encoded protein shares amino acid identity with ribosomal protein L23 from <i>Halobacterium marismortui</i>; however, its official symbol is RPL17. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring downstream C18orf32 (chromosome 18 open reading frame 32) gene. [provided by RefSeq, Dec 2010]</p>
<b>shRNA Design:</b>	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a>.</p>
<b>Performance Guaranteed:</b>	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>