

## Product datasheet for **TL309756**

### RPL13A Human shRNA Plasmid Kit (Locus ID 23521)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | RPL13A Human shRNA Plasmid Kit (Locus ID 23521)  |
| Locus ID:                 | 23521  |
| Synonyms:                 | L13A; TSTA1  |
| Vector:                   | pGFP-C-shLenti (TR30023)   |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)   |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Lentiviral plasmids  |
| Components:               | RPL13A - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 23521).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.  |
| RefSeq:                   | <a href="#">NM_001270491</a> , <a href="#">NM_012423</a> , <a href="#">NR_073024</a> , <a href="#">NM_012423.1</a> , <a href="#">NM_012423.2</a> , <a href="#">NM_012423.3</a> ,<br><a href="#">NM_001270491.1</a> , <a href="#">BC070223</a> , <a href="#">BC070223.1</a> , <a href="#">BC000514</a> , <a href="#">BC000847</a> , <a href="#">BC001675</a> , <a href="#">BC001836</a> ,<br><a href="#">BC004900</a> , <a href="#">BC013230</a> , <a href="#">BC032107</a> , <a href="#">BC062537</a> , <a href="#">BC065236</a> , <a href="#">BC071929</a> , <a href="#">BC105605</a> , <a href="#">NM_012423.4</a>   |
| UniProt ID:               | <a href="#">P40429</a>   |
| Summary:                  | Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a member of the L13P family of ribosomal proteins that is a component of the 60S subunit. The encoded protein also plays a role in the repression of inflammatory genes as a component of the IFN-gamma-activated inhibitor of translation (GAIT) complex. This gene is co-transcribed with the small nucleolar RNA genes U32, U33, U34, and U35, which are located in the second, fourth, fifth, and sixth introns, respectively. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed throughout the genome. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).