

Product datasheet for **TL309771V**

ROR gamma (RORC) Human shRNA Lentiviral Particle (Locus ID 6097)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	ROR gamma (RORC) Human shRNA Lentiviral Particle (Locus ID 6097)
Locus ID:	6097
Synonyms:	IMD42; NR1F3; RORG; RZR-GAMMA; RZRG; TOR
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	RORC - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC031554 , NM_001001523 , NM_005060 , NM_005060.1 , NM_005060.2 , NM_005060.3 , NM_001001523.1 , BC031554.1 , BC110571 , NM_001001523.2 , NM_005060.4
UniProt ID:	P51449
Summary:	The protein encoded by this gene is a DNA-binding transcription factor and is a member of the NR1 subfamily of nuclear hormone receptors. The specific functions of this protein are not known; however, studies of a similar gene in mice have shown that this gene may be essential for lymphoid organogenesis and may play an important regulatory role in thymopoiesis. In addition, studies in mice suggest that the protein encoded by this gene may inhibit the expression of Fas ligand and IL2. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).