

Product datasheet for **TL309831**

RGS6 Human shRNA Plasmid Kit (Locus ID 9628)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RGS6 Human shRNA Plasmid Kit (Locus ID 9628)
Locus ID:	9628
Synonyms:	GAP; HA117; S914
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RGS6 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9628). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<u>NM_001204416</u> , <u>NM_001204417</u> , <u>NM_001204418</u> , <u>NM_001204419</u> , <u>NM_001204420</u> , <u>NM_001204421</u> , <u>NM_001204422</u> , <u>NM_001204423</u> , <u>NM_001204424</u> , <u>NM_004296</u> , <u>NR_135235</u> , <u>NM_004296.1</u> , <u>NM_004296.2</u> , <u>NM_004296.3</u> , <u>NM_004296.4</u> , <u>NM_004296.5</u> , <u>NM_004296.6</u> , <u>NM_001204422.1</u> , <u>NM_001204422.2</u> , <u>NM_001204420.1</u> , <u>NM_001204420.2</u> , <u>NM_001204423.1</u> , <u>NM_001204421.1</u> , <u>NM_001204421.2</u> , <u>NM_001204419.1</u> , <u>NM_001204419.2</u> , <u>NM_001204418.1</u> , <u>NM_001204418.2</u> , <u>NM_001204417.1</u> , <u>NM_001204417.2</u> , <u>NM_001204416.1</u> , <u>NM_001204416.2</u> , <u>NM_001204424.1</u> , <u>BC104877</u> , <u>BC104877.1</u> , <u>BC113572</u> , <u>BC143699</u> , <u>NM_001370270</u> , <u>NM_001370271</u> , <u>NM_001370272</u> , <u>NM_001370273</u> , <u>NM_001370274</u> , <u>NM_001370276</u> , <u>NM_001370278</u> , <u>NM_001370281</u> , <u>NM_001370284</u> , <u>NM_001370290</u> , <u>NM_001370275</u> , <u>NM_001370277</u> , <u>NM_001370279</u> , <u>NM_001370280</u> , <u>NM_001370282</u> , <u>NM_001370283</u> , <u>NM_001370286</u> , <u>NM_001370287</u> , <u>NM_001370288</u> , <u>NM_001370289</u> , <u>NM_001370291</u> , <u>NM_001370292</u> , <u>NM_001370293</u> , <u>NM_001370294</u> , <u>NM_001204416.3</u> , <u>NM_001204423.2</u> , <u>NM_001204418.3</u> , <u>NM_001204417.3</u> , <u>NM_001204422.3</u> , <u>NM_001204424.2</u> , <u>NM_001204421.3</u> , <u>NM_004296.7</u> , <u>NM_001204420.3</u> , <u>NM_001204419.3</u>
UniProt ID:	<u>P49758</u>



[View online »](#)

Summary:	<p>This gene encodes a member of the RGS (regulator of G protein signaling) family of proteins, which are defined by the presence of a RGS domain that confers the GTPase-activating activity of these proteins toward certain G alpha subunits. This protein also belongs to a subfamily of RGS proteins characterized by the presence of DEP and GGL domains, the latter a G beta 5-interacting domain. The RGS proteins negatively regulate G protein signaling, and may modulate neuronal, cardiovascular, lymphocytic activities, and cancer risk. Many alternatively spliced transcript variants encoding different isoforms with long or short N-terminal domains, complete or incomplete GGL domains, and distinct C-terminal domains, have been described for this gene, however, the full-length nature of some of these variants is not known.[provided by RefSeq, Mar 2011]</p>
shRNA Design:	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.</p>
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>