

## Product datasheet for **TL309884V**

### Radixin (RDX) Human shRNA Lentiviral Particle (Locus ID 5962)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Radixin (RDX) Human shRNA Lentiviral Particle (Locus ID 5962)
Locus ID:	5962
Synonyms:	DFNB24
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	RDX - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001260492</a> , <a href="#">NM_001260493</a> , <a href="#">NM_001260494</a> , <a href="#">NM_001260495</a> , <a href="#">NM_001260496</a> , <a href="#">NM_002906</a> , <a href="#">NM_002906.1</a> , <a href="#">NM_002906.2</a> , <a href="#">NM_002906.3</a> , <a href="#">NM_001260496.1</a> , <a href="#">NM_001260495.1</a> , <a href="#">NM_001260494.1</a> , <a href="#">NM_001260492.1</a> , <a href="#">NM_001260493.1</a> , <a href="#">BC047109</a> , <a href="#">BC002626</a> , <a href="#">BC020751</a> , <a href="#">BC029467</a> , <a href="#">NM_001260492.2</a> , <a href="#">NM_002906.4</a>
UniProt ID:	<a href="#">P35241</a>
Summary:	Radixin is a cytoskeletal protein that may be important in linking actin to the plasma membrane. It is highly similar in sequence to both ezrin and moesin. The radixin gene has been localized by fluorescence in situ hybridization to 11q23. A truncated version representing a pseudogene (RDXP2) was assigned to Xp21.3. Another pseudogene that seemed to lack introns (RDXP1) was mapped to 11p by Southern and PCR analyses. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).