

Product datasheet for **TL309886V**

Retinol dehydrogenase 16 (RDH16) Human shRNA Lentiviral Particle (Locus ID 8608)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Retinol dehydrogenase 16 (RDH16) Human shRNA Lentiviral Particle (Locus ID 8608)
Locus ID:	8608
Synonyms:	hRDH-E; RODH-4; SDR9C8
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	RDH16 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001320108 , NM_003708 , NM_003708.1 , NM_003708.2 , NM_003708.3 , NM_003708.4 , BC160081
UniProt ID:	O75452
Summary:	Oxidoreductase with a preference for NAD. Oxidizes all-trans-retinol, 9-cis-retinol, 11-cis-retinol and 13-cis-retinol to the corresponding aldehydes (PubMed:10329026, PubMed:12534290, PubMed:9677409). Has higher activity towards CRBP-bound retinol than with free retinol (PubMed:12534290). Oxidizes also 3-alpha-hydroxysteroids. Oxidizes androstenediol and androsterone to dihydrotestosterone and androstenedione. Can also catalyze the reverse reaction (PubMed:10329026, PubMed:9677409, PubMed:29541409). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).