

## Product datasheet for **TL309923**

### **RbAp48 (RBBP4) Human shRNA Plasmid Kit (Locus ID 5928)**

#### **Product data:**

|                                  |  |
|----------------------------------|--|
| <b>Product Type:</b>             | shRNA Plasmids   |
| <b>Product Name:</b>             | RbAp48 (RBBP4) Human shRNA Plasmid Kit (Locus ID 5928)   |
| <b>Locus ID:</b>                 | 5928   |
| <b>Synonyms:</b>                 | lin-53; NURF55; RBAP48   |
| <b>Vector:</b>                   | pGFP-C-shLenti (TR30023)   |
| <b>E. coli Selection:</b>        | Chloramphenicol (34 ug/ml)   |
| <b>Mammalian Cell Selection:</b> | Puromycin  |
| <b>Format:</b>                   | Lentiviral plasmids  |
| <b>Components:</b>               | RBBP4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 5928).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.  |
| <b>RefSeq:</b>                   | <a href="#">NM_001135255</a> , <a href="#">NM_001135256</a> , <a href="#">NM_005610</a> , <a href="#">NM_005610.1</a> , <a href="#">NM_005610.2</a> , <a href="#">NM_001135256.1</a> , <a href="#">NM_001135255.1</a> , <a href="#">BC053904</a> , <a href="#">BC053904.1</a> , <a href="#">BC075836</a> , <a href="#">BC003092</a> , <a href="#">BC015123</a> , <a href="#">NM_001135255.2</a>  |
| <b>UniProt ID:</b>               | <a href="#">Q09028</a>   |
| <b>Summary:</b>                  | This gene encodes a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. It is present in protein complexes involved in histone acetylation and chromatin assembly. It is part of the Mi-2 complex which has been implicated in chromatin remodeling and transcriptional repression associated with histone deacetylation. This encoded protein is also part of co-repressor complexes, which is an integral component of transcriptional silencing. It is found among several cellular proteins that bind directly to retinoblastoma protein to regulate cell proliferation. This protein also seems to be involved in transcriptional repression of E2F-responsive genes. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2008] |
| <b>shRNA Design:</b>             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).