

Product datasheet for **TL309968V**

RAD21 Human shRNA Lentiviral Particle (Locus ID 5885)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	RAD21 Human shRNA Lentiviral Particle (Locus ID 5885)
Locus ID:	5885
Synonyms:	CDLS4; hHR21; HR21; HRAD21; MCD1; MGS; NXP1; SCC1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	RAD21 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, $>10^7$ TU/ml.
RefSeq:	NM_006265 , NM_006265.1 , NM_006265.2 , BC050381 , BC050381.1 , BC001229 , NM_006265.3
UniProt ID:	Q60216
Summary:	The protein encoded by this gene is highly similar to the gene product of <i>Schizosaccharomyces pombe rad21</i> , a gene involved in the repair of DNA double-strand breaks, as well as in chromatid cohesion during mitosis. This protein is a nuclear phospho-protein, which becomes hyperphosphorylated in cell cycle M phase. The highly regulated association of this protein with mitotic chromatin specifically at the centromere region suggests its role in sister chromatid cohesion in mitotic cells. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).