

Product datasheet for **TL309976**

RABGAP1L Human shRNA Plasmid Kit (Locus ID 9910)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RABGAP1L Human shRNA Plasmid Kit (Locus ID 9910)
Locus ID:	9910
Synonyms:	HHL; TBC1D18
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RABGAP1L - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9910). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001035230 , NM_001243763 , NM_001243764 , NM_001243765 , NM_014857 , NM_001330989 , NM_014857.1 , NM_014857.2 , NM_014857.3 , NM_014857.4 , NM_001035230.1 , NM_001035230.2 , NM_001243763.1 , NM_001243764.1 , NM_001243765.1 , BC012094 , BC018630 , BC041888 , BC062541 , NM_001366445 , NM_001366447 , NM_001366448 , NM_001366451 , NM_001366452 , NM_001366453 , NM_001366457 , NR_158982 , NM_001366446 , NM_001366449 , NM_001366450 , NM_001366454 , NM_001366455 , NM_001366456 , NR_158981 , NM_001243765.2 , NM_001243764.2
UniProt ID:	B7ZAP0
Summary:	GTP-hydrolysis activating protein (GAP) for small GTPase RAB22A, converting active RAB22A-GTP to the inactive form RAB22A-GDP (PubMed:16923123). Plays a role in endocytosis and intracellular protein transport. Recruited by ANK2 to phosphatidylinositol 3-phosphate (PI3P)-positive early endosomes, where it inactivates RAB22A, and promotes polarized trafficking to the leading edge of the migrating cells. Part of the ANK2/RABGAP1L complex which is required for the polarized recycling of fibronectin receptor ITGA5 ITGB1 to the plasma membrane that enables continuous directional cell migration (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).