

## Product datasheet for **TL310167**

### PRMT5 Human shRNA Plasmid Kit (Locus ID 10419)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	PRMT5 Human shRNA Plasmid Kit (Locus ID 10419)
Locus ID:	10419
Synonyms:	HRMT1L5; HSL7; IBP72; JBP1; SKB1; SKB1Hs
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PRMT5 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10419). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001039619</a> , <a href="#">NM_001282953</a> , <a href="#">NM_001282954</a> , <a href="#">NM_001282955</a> , <a href="#">NM_001282956</a> , <a href="#">NM_006109</a> , <a href="#">NM_006109.1</a> , <a href="#">NM_006109.2</a> , <a href="#">NM_006109.3</a> , <a href="#">NM_006109.4</a> , <a href="#">NM_001039619.1</a> , <a href="#">NM_001039619.2</a> , <a href="#">NM_001282956.1</a> , <a href="#">NM_001282954.1</a> , <a href="#">NM_001282953.1</a> , <a href="#">NM_001282955.1</a> , <a href="#">BC025979</a> , <a href="#">BC025979.2</a> , <a href="#">BC005820</a> , <a href="#">NM_001282956.2</a> , <a href="#">NM_001039619.3</a> , <a href="#">NM_001282954.2</a> , <a href="#">NM_001282953.2</a> , <a href="#">NM_006109.5</a> , <a href="#">NM_001282955.2</a>
UniProt ID:	<a href="#">Q14744</a>
Summary:	This gene encodes an enzyme that belongs to the methyltransferase family. The encoded protein catalyzes the transfer of methyl groups to the amino acid arginine, in target proteins that include histones, transcriptional elongation factors and the tumor suppressor p53. This gene plays a role in several cellular processes, including transcriptional regulation, and the assembly of small nuclear ribonucleoproteins. A pseudogene of this gene has been defined on chromosome 4. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).