

Product datasheet for **TL310215**

PP5 (PPP5C) Human shRNA Plasmid Kit (Locus ID 5536)

Product data:

Product Type:	shRNA Plasmids
Product Name:	PP5 (PPP5C) Human shRNA Plasmid Kit (Locus ID 5536)
Locus ID:	5536
Synonyms:	PP5; PPP5; PPT
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PPP5C - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 5536). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001204284 , NM_006247 , NM_006247.1 , NM_006247.2 , NM_006247.3 , NM_001204284.1 , BC001970 , BC001970.1 , BC000750 , BC001831 , NM_006247.4 , NM_001204284.2
UniProt ID:	P53041
Summary:	This gene encodes a serine/threonine phosphatase which is a member of the protein phosphatase catalytic subunit family. Proteins in this family participate in pathways regulated by reversible phosphorylation at serine and threonine residues; many of these pathways are involved in the regulation of cell growth and differentiation. The product of this gene has been shown to participate in signaling pathways in response to hormones or cellular stress, and elevated levels of this protein may be associated with breast cancer development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL310215A virus into HEK293 cells. TL310215A virus was prepared using lenti-shRNA TL310215A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL310215B virus into HEK293 cells. TL310215B virus was prepared using lenti-shRNA TL310215B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL310215C] virus into HEK293 cells. [TL310215C] virus was prepared using lenti-shRNA [TL310215C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL310215D] virus into HEK293 cells. [TL310215D] virus was prepared using lenti-shRNA [TL310215D] and [TR30037] packaging kit.