

## Product datasheet for **TL310226**

### PPP2R1B Human shRNA Plasmid Kit (Locus ID 5519)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	PPP2R1B Human shRNA Plasmid Kit (Locus ID 5519)
Locus ID:	5519
Synonyms:	PP2A-Abeta; PR65B
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PPP2R1B - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 5519). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001177562</a> , <a href="#">NM_001177563</a> , <a href="#">NM_002716</a> , <a href="#">NM_181699</a> , <a href="#">NM_181700</a> , <a href="#">NM_181699.1</a> , <a href="#">NM_181699.2</a> , <a href="#">NM_002716.1</a> , <a href="#">NM_002716.2</a> , <a href="#">NM_002716.3</a> , <a href="#">NM_002716.4</a> , <a href="#">NM_001177563.1</a> , <a href="#">NM_001177562.1</a> , <a href="#">NM_181700.1</a> , <a href="#">BC027596</a> , <a href="#">BC027596.2</a> , <a href="#">NM_002716.5</a> , <a href="#">NM_181700.2</a> , <a href="#">NM_001177563.2</a> , <a href="#">NM_001177562.2</a> , <a href="#">NM_181699.3</a>
UniProt ID:	<a href="#">P30154</a>
Summary:	This gene encodes a constant regulatory subunit of protein phosphatase 2. Protein phosphatase 2 is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The constant regulatory subunit A serves as a scaffolding molecule to coordinate the assembly of the catalytic subunit and a variable regulatory B subunit. This gene encodes a beta isoform of the constant regulatory subunit A. Mutations in this gene have been associated with some lung and colon cancers. Alternatively spliced transcript variants have been described. [provided by RefSeq, Apr 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).