

Product datasheet for **TL310381**

PLAGL1 Human shRNA Plasmid Kit (Locus ID 5325)

Product data:

Product Type:	shRNA Plasmids
Product Name:	PLAGL1 Human shRNA Plasmid Kit (Locus ID 5325)
Locus ID:	5325
Synonyms:	LOT1; ZAC; ZAC1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PLAGL1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 5325). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001080951 , NM_001080952 , NM_001080953 , NM_001080954 , NM_001080955 , NM_001080956 , NM_001289037 , NM_001289038 , NM_001289039 , NM_001289040 , NM_001289041 , NM_001289042 , NM_001289043 , NM_001289044 , NM_001289045 , NM_001289046 , NM_001289047 , NM_001289048 , NM_001289049 , NM_001317156 , NM_001317157 , NM_001317158 , NM_001317159 , NM_001317160 , NM_001317161 , NM_001317162 , NM_002656 , NM_006718 , NM_001080954.1 , NM_001080954.2 , NM_001080955.1 , NM_001080955.2 , NM_001080952.1 , NM_001080952.2 , NM_006718.1 , NM_006718.2 , NM_006718.3 , NM_006718.4 , NM_001080951.1 , NM_001080951.2 , NM_001080953.1 , NM_001080953.2 , NM_001080956.1 , NM_001080956.2 , NM_001289037.1 , NM_001289038.1 , NM_001289039.1 , NM_001289040.1 , NM_001289041.1 , NM_001289042.1 , NM_001289043.1 , NM_001289044.1 , NM_001289045.1 , NM_001289046.1 , NM_001289047.1 , NM_001289048.1 , NM_001289049.1 , NM_002656.2 , NM_002656.3 , BC109086 , BC074814 , BC109085 , NM_006718.5 , NM_001289039.2 , NM_001289037.2 , NM_001080954.3 , NM_001289044.2 , NM_001080953.3 , NM_001080955.3 , NM_001289040.2 , NM_001289049.2 , NM_001289041.2 , NM_001080951.3 , NM_001289048.2 , NM_001080956.3 , NM_001289046.2 , NM_001289042.2 , NM_001289045.2 , NM_001080952.3
UniProt ID:	Q9UM63



[View online »](#)

Summary:	<p>This gene encodes a C2H2 zinc finger protein that functions as a suppressor of cell growth. This gene is often deleted or methylated and silenced in cancer cells. In addition, overexpression of this gene during fetal development is thought to be the causal factor for transient neonatal diabetes mellitus (TNDM). Alternative splicing and the use of alternative promoters results in multiple transcript variants encoding two different protein isoforms. The P1 downstream promoter of this gene is imprinted, with preferential expression from the paternal allele in many tissues. [provided by RefSeq, Nov 2015]</p>
shRNA Design:	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.</p>
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>