

## Product datasheet for **TL310484V**

### PFKFB3 Human shRNA Lentiviral Particle (Locus ID 5209)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	PFKFB3 Human shRNA Lentiviral Particle (Locus ID 5209)
Locus ID:	5209
Synonyms:	iPFK-2; IPFK2; PFK2
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	PFKFB3 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001145443</a> , <a href="#">NM_001282630</a> , <a href="#">NM_001314063</a> , <a href="#">NM_001323016</a> , <a href="#">NM_001323017</a> , <a href="#">NM_004566</a> , <a href="#">NR_136554</a> , <a href="#">NM_004566.1</a> , <a href="#">NM_004566.2</a> , <a href="#">NM_004566.3</a> , <a href="#">NM_001145443.1</a> , <a href="#">NM_001145443.2</a> , <a href="#">NM_001282630.1</a> , <a href="#">NM_001282630.2</a> , <a href="#">BC040482</a> , <a href="#">BC040482.1</a> , <a href="#">BC037258</a> , <a href="#">BC042656</a> , <a href="#">NM_001363545</a> , <a href="#">NM_004566.4</a> , <a href="#">NM_001145443.3</a>
UniProt ID:	<a href="#">Q16875</a>
Summary:	The protein encoded by this gene belongs to a family of bifunctional proteins that are involved in both the synthesis and degradation of fructose-2,6-bisphosphate, a regulatory molecule that controls glycolysis in eukaryotes. The encoded protein has a 6-phosphofructo-2-kinase activity that catalyzes the synthesis of fructose-2,6-bisphosphate (F2,6BP), and a fructose-2,6-bisphosphatase activity that catalyzes the degradation of F2,6BP. This protein is required for cell cycle progression and prevention of apoptosis. It functions as a regulator of cyclin-dependent kinase 1, linking glucose metabolism to cell proliferation and survival in tumor cells. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL310484A virus into HEK293 cells. TL310484A virus was prepared using lenti-shRNA TL310484A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL310484B virus into HEK293 cells. TL310484B virus was prepared using lenti-shRNA TL310484B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL310484C] virus into HEK293 cells. [TL310484C] virus was prepared using lenti-shRNA [TL310484C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL310484D] virus into HEK293 cells. [TL310484D] virus was prepared using lenti-shRNA [TL310484D] and [TR30037] packaging kit.