

## Product datasheet for **TL310599V**

### **PAX2 Human shRNA Lentiviral Particle (Locus ID 5076)**

#### **Product data:**

<b>Product Type:</b>	shRNA Lentiviral Particles
<b>Product Name:</b>	PAX2 Human shRNA Lentiviral Particle (Locus ID 5076)
<b>Locus ID:</b>	5076
<b>Synonyms:</b>	FSGS7; PAPRS
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>Format:</b>	Lentiviral particles
<b>Components:</b>	PAX2 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
<b>RefSeq:</b>	<a href="#">NM_000278</a> , <a href="#">NM_001304569</a> , <a href="#">NM_003987</a> , <a href="#">NM_003988</a> , <a href="#">NM_003989</a> , <a href="#">NM_003990</a> , <a href="#">NM_003990.1</a> , <a href="#">NM_003990.2</a> , <a href="#">NM_003990.3</a> , <a href="#">NM_003990.4</a> , <a href="#">NM_003987.1</a> , <a href="#">NM_003987.3</a> , <a href="#">NM_000278.1</a> , <a href="#">NM_000278.2</a> , <a href="#">NM_000278.3</a> , <a href="#">NM_003988.2</a> , <a href="#">NM_003988.3</a> , <a href="#">NM_003988.4</a> , <a href="#">NM_003989.1</a> , <a href="#">NM_003989.2</a> , <a href="#">NM_003989.3</a> , <a href="#">NM_003989.4</a> , <a href="#">BC141452</a> , <a href="#">BC148710</a> , <a href="#">BM671839</a> , <a href="#">NM_003988.5</a> , <a href="#">NM_000278.5</a>
<b>UniProt ID:</b>	<a href="#">Q02962</a>
<b>Summary:</b>	PAX2 encodes paired box gene 2, one of many human homologues of the Drosophila melanogaster gene prd. The central feature of this transcription factor gene family is the conserved DNA-binding paired box domain. PAX2 is believed to be a target of transcriptional suppression by the tumor suppressor gene WT1. Mutations within PAX2 have been shown to result in optic nerve colobomas and renal hypoplasia. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2014]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL310599A virus into HEK293 cells. TL310599A virus was prepared using lenti-shRNA TL310599A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL310599B virus into HEK293 cells. TL310599B virus was prepared using lenti-shRNA TL310599B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL310599C] virus into HEK293 cells. [TL310599C] virus was prepared using lenti-shRNA [TL310599C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL310599D] virus into HEK293 cells. [TL310599D] virus was prepared using lenti-shRNA [TL310599D] and [TR30037] packaging kit.