

## Product datasheet for **TL310610**

### **PARG Human shRNA Plasmid Kit (Locus ID 8505)**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | PARG Human shRNA Plasmid Kit (Locus ID 8505)   |
| Locus ID:                 | 8505   |
| Synonyms:                 | PARG99   |
| Vector:                   | pGFP-C-shLenti (TR30023)   |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)   |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Lentiviral plasmids  |
| Components:               | PARG - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 8505).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.   |
| RefSeq:                   | <a href="#">NM_001303486</a> , <a href="#">NM_001303487</a> , <a href="#">NM_001303489</a> , <a href="#">NM_001324381</a> , <a href="#">NM_003631</a> , <a href="#">NR_130168</a> , <a href="#">NR_130169</a> , <a href="#">NR_136752</a> , <a href="#">NR_136753</a> , <a href="#">NR_136754</a> , <a href="#">NR_136755</a> , <a href="#">NR_136756</a> , <a href="#">NM_003631.1</a> , <a href="#">NM_003631.2</a> , <a href="#">NM_003631.3</a> , <a href="#">BC050560</a> , <a href="#">BC012005</a> , <a href="#">BC052966</a> , <a href="#">NM_003631.5</a> |
| UniProt ID:               | <a href="#">Q86W56</a>   |
| Summary:                  | Poly(ADP-ribose) glycohydrolase (PARG) is the major enzyme responsible for the catabolism of poly(ADP-ribose), a reversible covalent-modifier of chromosomal proteins. The protein is found in many tissues and may be subject to proteolysis generating smaller, active products. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2015]   |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL310610A virus into HEK293 cells. TL310610A virus was prepared using lenti-shRNA TL310610A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL310610B virus into HEK293 cells. TL310610B virus was prepared using lenti-shRNA TL310610B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL310610C] virus into HEK293 cells. [TL310610C] virus was prepared using lenti-shRNA [TL310610C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL310610D] virus into HEK293 cells. [TL310610D] virus was prepared using lenti-shRNA [TL310610D] and [TR30037] packaging kit.