

## Product datasheet for **TL311096**

### NRG2 Human shRNA Plasmid Kit (Locus ID 9542)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	NRG2 Human shRNA Plasmid Kit (Locus ID 9542)
Locus ID:	9542
Synonyms:	DON1; HRG2; NTAK
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	NRG2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9542). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001184935</a> , <a href="#">NM_004883</a> , <a href="#">NM_013981</a> , <a href="#">NM_013982</a> , <a href="#">NM_013983</a> , <a href="#">NM_013984</a> , <a href="#">NM_013985</a> , <a href="#">NM_013981.1</a> , <a href="#">NM_013981.2</a> , <a href="#">NM_013981.3</a> , <a href="#">NM_013982.1</a> , <a href="#">NM_013982.2</a> , <a href="#">NM_004883.1</a> , <a href="#">NM_004883.2</a> , <a href="#">NM_013983.1</a> , <a href="#">NM_013983.2</a> , <a href="#">NM_001184935.1</a> , <a href="#">BC103983</a> , <a href="#">BC103984</a> , <a href="#">BC114333</a> , <a href="#">BC114334</a> , <a href="#">BC166615</a> , <a href="#">NM_004883.3</a> , <a href="#">NM_013982.3</a> , <a href="#">NM_013983.3</a> , <a href="#">NM_013981.4</a>
UniProt ID:	<a href="#">O14511</a>
Summary:	This gene encodes a novel member of the neuregulin family of growth and differentiation factors. Through interaction with the ERBB family of receptors, this protein induces the growth and differentiation of epithelial, neuronal, glial, and other types of cells. The gene consists of 12 exons and the genomic structure is similar to that of neuregulin 1, another member of the neuregulin family of ligands. The products of these genes mediate distinct biological processes by acting at different sites in tissues and eliciting different biological responses in cells. This gene is located close to the region for demyelinating Charcot-Marie-Tooth disease locus, but is not responsible for this disease. Alternative transcript variants encoding distinct isoforms have been described. [provided by RefSeq, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).