

## Product datasheet for **TL311182**

### NFYA Human shRNA Plasmid Kit (Locus ID 4800)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	NFYA Human shRNA Plasmid Kit (Locus ID 4800)
Locus ID:	4800
Synonyms:	CBF-A; CBF-B; HAP2; NF-YA
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	NFYA - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 4800). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_002505</a> , <a href="#">NM_021705</a> , <a href="#">NM_021705.1</a> , <a href="#">NM_021705.2</a> , <a href="#">NM_021705.3</a> , <a href="#">NM_002505.1</a> , <a href="#">NM_002505.2</a> , <a href="#">NM_002505.3</a> , <a href="#">NM_002505.4</a> , <a href="#">BC039244</a> , <a href="#">BC039244.1</a> , <a href="#">BC026923</a> , <a href="#">BM561405</a> , <a href="#">NM_002505.5</a> , <a href="#">NM_021705.4</a>
UniProt ID:	<a href="#">P23511</a>
Summary:	The protein encoded by this gene is one subunit of a trimeric complex, forming a highly conserved transcription factor that binds to CCAAT motifs in the promoter regions in a variety of genes. Subunit A associates with a tight dimer composed of the B and C subunits, resulting in a trimer that binds to DNA with high specificity and affinity. The sequence specific interactions of the complex are made by the A subunit, suggesting a role as the regulatory subunit. In addition, there is evidence of post-transcriptional regulation in this gene product, either by protein degradation or control of translation. Further regulation is represented by alternative splicing in the glutamine-rich activation domain, with clear tissue-specific preferences for the two isoforms. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).