

Product datasheet for TL311194

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Nrf2 (NFE2L2) Human shRNA Plasmid Kit (Locus ID 4780)

Product data:

Product Type: shRNA Plasmids

Product Name: Nrf2 (NFE2L2) Human shRNA Plasmid Kit (Locus ID 4780)

Locus ID: 4780

Synonyms: HEBP1; IMDDHH; Nrf-2; NRF2

Vector: pGFP-C-shLenti (TR30023)

E. coli Selection: Chloramphenicol (34 ug/ml)

Mammalian Cell

Selection:

Puromycin

Format: Lentiviral plasmids

Components: NFE2L2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 4780).

5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.

RefSeq: NM 001145412, NM 001145413, NM 001313900, NM 001313901, NM 001313902,

NM 001313903, NM 001313904, NM 006164, NM 006164.1, NM 006164.2, NM 006164.3, NM 006164.4, NM 001145412.1, NM 001145412.2, NM 001145412.3, NM 001145413.1,

NM 001145413.2, NM 001145413.3, BC011558, BC011558.1, NM 006164.5

UniProt ID: Q16236

Summary: This gene encodes a transcription factor which is a member of a small family of basic leucine

zipper (bZIP) proteins. The encoded transcription factor regulates genes which contain antioxidant response elements (ARE) in their promoters; many of these genes encode proteins involved in response to injury and inflammation which includes the production of free radicals. Multiple transcript variants encoding different isoforms have been characterized

for this gene. [provided by RefSeq, Sep 2015]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u>.






Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:




GFP signal was observed under microscope at 48 hours after transduction of TL311194A virus into HEK293 cells. TL311194A virus was prepared using lenti-shRNA TL311194A and [TR30037] packaging kit.






GFP signal was observed under microscope at 48 hours after transduction of TL311194B virus into HEK293 cells. TL311194B virus was prepared using lenti-shRNA TL311194B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL311194C] virus into HEK293 cells. [TL311194C] virus was prepared using lenti-shRNA [TL311194C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL311194D] virus into HEK293 cells. [TL311194D] virus was prepared using lenti-shRNA [TL311194D] and [TR30037] packaging kit.