

## Product datasheet for **TL311199V**

### Neurofibromin (NF1) Human shRNA Lentiviral Particle (Locus ID 4763)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Neurofibromin (NF1) Human shRNA Lentiviral Particle (Locus ID 4763)
Locus ID:	4763
Synonyms:	NFNS; VRNF; WSS
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	NF1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<u><a href="#">NM_000267</a></u> , <u><a href="#">NM_001042492</a></u> , <u><a href="#">NM_001128147</a></u> , <u><a href="#">NM_001042492.1</a></u> , <u><a href="#">NM_001042492.2</a></u> , <u><a href="#">NM_001128147.1</a></u> , <u><a href="#">NM_001128147.2</a></u> , <u><a href="#">NM_000267.1</a></u> , <u><a href="#">NM_000267.3</a></u> , <u><a href="#">BC144643</a></u> , <u><a href="#">NM_001128147.3</a></u> , <u><a href="#">NM_001042492.3</a></u>
UniProt ID:	<u><a href="#">P21359</a></u>
Summary:	This gene product appears to function as a negative regulator of the ras signal transduction pathway. Mutations in this gene have been linked to neurofibromatosis type 1, juvenile myelomonocytic leukemia and Watson syndrome. The mRNA for this gene is subject to RNA editing (CGA>UGA->Arg1306Term) resulting in premature translation termination. Alternatively spliced transcript variants encoding different isoforms have also been described for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).