

Product datasheet for **TL311300**

MYO1A Human shRNA Plasmid Kit (Locus ID 4640)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MYO1A Human shRNA Plasmid Kit (Locus ID 4640)
Locus ID:	4640
Synonyms:	BBMI; DFNA48; MIHC; MYHL
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MYO1A - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 4640). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001256041 , NM_005379 , NM_005379.1 , NM_005379.2 , NM_005379.3 , NM_001256041.1 , BC059387 , BC059387.1 , NM_001256041.2 , NM_005379.4
UniProt ID:	Q9UBC5
Summary:	This gene encodes a member of the myosin superfamily. The protein represents an unconventional myosin; it should not be confused with the conventional skeletal muscle myosin-1 (MYH1). Unconventional myosins contain the basic domains characteristic of conventional myosins and are further distinguished from class members by their tail domains. They function as actin-based molecular motors. Mutations in this gene have been associated with autosomal dominant deafness. Alternatively spliced variants have been found for this gene. [provided by RefSeq, Dec 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).