

Product datasheet for **TL311343**

MTUS1 Human shRNA Plasmid Kit (Locus ID 57509)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MTUS1 Human shRNA Plasmid Kit (Locus ID 57509)
Locus ID:	57509
Synonyms:	ATBP; ATIP; ATIP3; ICIS; MP44; MTSG1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MTUS1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 57509). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001001924 , NM_001001925 , NM_001001927 , NM_001001931 , NM_001166393 , NM_001330470 , NM_020749 , NM_001001924.1 , NM_020749.1 , NM_020749.2 , NM_020749.3 , NM_020749.4 , NM_001001925.1 , NM_001001925.2 , NM_001001931.1 , NM_001001931.2 , NM_001166393.1 , NM_001001927.1 , BC007328 , BC017740 , BC033842 , BC098381 , BC117564 , BC136320 , BC140067 , BC140068 , BC142971 , BC171765 , NM_001363057 , NM_001363059 , NM_001363061 , NM_001363062 , NM_001363058 , NM_001363060 , NR_156437
UniProt ID:	Q9ULD2
Summary:	This gene encodes a protein which contains a C-terminal domain able to interact with the angiotension II (AT2) receptor and a large coiled-coil region allowing dimerization. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. One of the transcript variants has been shown to encode a mitochondrial protein that acts as a tumor suppressor and participates in AT2 signaling pathways. Other variants may encode nuclear or transmembrane proteins but it has not been determined whether they also participate in AT2 signaling pathways. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).