

## Product datasheet for **TL311380**

### MSH6 Human shRNA Plasmid Kit (Locus ID 2956)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	MSH6 Human shRNA Plasmid Kit (Locus ID 2956)
Locus ID:	2956
Synonyms:	GTBP; GTMBP; HNPCC5; HSAP; MMRCS3; p160
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MSH6 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 2956). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_000179</a> , <a href="#">NM_001281492</a> , <a href="#">NM_001281493</a> , <a href="#">NM_001281494</a> , <a href="#">NM_000179.1</a> , <a href="#">NM_000179.2</a> , <a href="#">NM_001281493.1</a> , <a href="#">NM_001281494.1</a> , <a href="#">NM_001281492.1</a> , <a href="#">BC004246</a> , <a href="#">BC004246.1</a> , <a href="#">BC071594</a> , <a href="#">BC104665</a> , <a href="#">NM_000179.3</a>
UniProt ID:	<a href="#">P52701</a>
Summary:	This gene encodes a member of the DNA mismatch repair MutS family. In E. coli, the MutS protein helps in the recognition of mismatched nucleotides prior to their repair. A highly conserved region of approximately 150 aa, called the Walker-A adenine nucleotide binding motif, exists in MutS homologs. The encoded protein heterodimerizes with MSH2 to form a mismatch recognition complex that functions as a bidirectional molecular switch that exchanges ADP and ATP as DNA mismatches are bound and dissociated. Mutations in this gene may be associated with hereditary nonpolyposis colon cancer, colorectal cancer, and endometrial cancer. Transcripts variants encoding different isoforms have been described. [provided by RefSeq, Jul 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).