

Product datasheet for **TL311397**

MRPL18 Human shRNA Plasmid Kit (Locus ID 29074)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MRPL18 Human shRNA Plasmid Kit (Locus ID 29074)
Locus ID:	29074
Synonyms:	HSPC071; L18mt; MRP-L18
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MRPL18 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 29074). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001318817 , NM_014161 , NR_134860 , NR_134861 , NR_134862 , NR_134863 , NR_134864 , NR_134865 , NR_134866 , NM_014161.1 , NM_014161.2 , NM_014161.3 , NM_014161.4 , BC001623 , BC001623.1 , BM924328 , NM_014161.5
UniProt ID:	Q9H0U6
Summary:	This nuclear gene encodes a protein component of the larger 39S subunit of mitochondrial ribosome. This protein may also aid in the import of nuclear-encoded 5S rRNA into mitochondria. Alternative splicing results in multiple transcript variants, most of which are not predicted to encode a protein. A pseudogene of this gene is found on chromosome 16. [provided by RefSeq, Jan 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).