

Product datasheet for **TL311481V**

Monoacylglycerol Lipase (MGLL) Human shRNA Lentiviral Particle (Locus ID 11343)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Monoacylglycerol Lipase (MGLL) Human shRNA Lentiviral Particle (Locus ID 11343)
Locus ID:	11343
Synonyms:	HU-K5; HUK5; MAGL; MGL
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	MGLL - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, $>10^7$ TU/ml.
RefSeq:	NM_001003794 , NM_001256585 , NM_007283 , NM_001003794.1 , NM_001003794.2 , NM_007283.1 , NM_007283.2 , NM_007283.3 , NM_007283.4 , NM_007283.5 , NM_007283.6 , NM_001256585.1 , BC006230 , BC006230.2 , BC000551 , BC047298 , BC073823 , BM051374 , BM669411
UniProt ID:	Q99685
Summary:	This gene encodes a serine hydrolase of the AB hydrolase superfamily that catalyzes the conversion of monoacylglycerides to free fatty acids and glycerol. The encoded protein plays a critical role in several physiological processes including pain and nociception through hydrolysis of the endocannabinoid 2-arachidonoylglycerol. Expression of this gene may play a role in cancer tumorigenesis and metastasis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Feb 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).