

Product datasheet for **TL311583V**

NFkB Inducing Kinase NIK (MAP3K14) Human shRNA Lentiviral Particle (Locus ID 9020)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	NFkB Inducing Kinase NIK (MAP3K14) Human shRNA Lentiviral Particle (Locus ID 9020)
Locus ID:	9020
Synonyms:	FTDCR1B; HS; HSNIK; NIK
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	MAP3K14 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC035576 , NM_003954 , NM_003954.1 , NM_003954.2 , NM_003954.3 , NM_003954.4 , BC035576.1 , BC022296 , NM_003954.5
UniProt ID:	Q99558
Summary:	This gene encodes mitogen-activated protein kinase kinase kinase 14, which is a serine/threonine protein-kinase. This kinase binds to TRAF2 and stimulates NF-kappaB activity. It shares sequence similarity with several other MAPKK kinases. It participates in an NF-kappaB-inducing signalling cascade common to receptors of the tumour-necrosis/nerve-growth factor (TNF/NGF) family and to the interleukin-1 type-I receptor. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:

GFP signal was observed under microscope at 48 hours after transduction of TL311583A virus into HEK293 cells. TL311583A virus was prepared using lenti-shRNA TL311583A and [TR30037] packaging kit.