

Product datasheet for **TL311614**

MAGEA2 Human shRNA Plasmid Kit (Locus ID 4101)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MAGEA2 Human shRNA Plasmid Kit (Locus ID 4101)
Locus ID:	4101
Synonyms:	CT1.2; MAGE2; MAGEA2A
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MAGEA2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 4101). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001282501 , NM_001282502 , NM_001282504 , NM_001282505 , NM_005361 , NM_175742 , NM_175743 , NM_175743.1 , NM_175743.2 , NM_175742.1 , NM_175742.2 , NM_005361.1 , NM_005361.2 , NM_005361.3 , NM_001282501.1 , NM_001282502.1 , NM_001282504.1 , NM_001282505.1 , BC063681 , BC108720 , BC013098 , BC109188 , BC112158 , BC112160 , BC146449 , BM781931
UniProt ID:	P43356
Summary:	This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. This gene has two identical copies at different loci. Alternatively spliced transcript variants encoding the same protein have been identified for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).