

## Product datasheet for **TL311737**

### LIFR Human shRNA Plasmid Kit (Locus ID 3977)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	LIFR Human shRNA Plasmid Kit (Locus ID 3977)
Locus ID:	3977
Synonyms:	CD118; LIF-R; SJS2; STWS; SWS
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	LIFR - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3977). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001127671</a> , <a href="#">NM_002310</a> , <a href="#">NM_002310.1</a> , <a href="#">NM_002310.2</a> , <a href="#">NM_002310.3</a> , <a href="#">NM_002310.4</a> , <a href="#">NM_002310.5</a> , <a href="#">NM_001127671.1</a> , <a href="#">BC141428</a> , <a href="#">BC153096</a> , <a href="#">NM_001364297</a> , <a href="#">NM_001364298</a> , <a href="#">NM_001127671.2</a> , <a href="#">NM_002310.6</a>
UniProt ID:	<a href="#">P42702</a>
Summary:	This gene encodes a protein that belongs to the type I cytokine receptor family. This protein combines with a high-affinity converter subunit, gp130, to form a receptor complex that mediates the action of the leukemia inhibitory factor, a polyfunctional cytokine that is involved in cellular differentiation, proliferation and survival in the adult and the embryo. Mutations in this gene cause Schwartz-Jampel syndrome type 2, a disease belonging to the group of the bent-bone dysplasias. A translocation that involves the promoter of this gene, t(5;8)(p13;q12) with the pleiomorphic adenoma gene 1, is associated with salivary gland pleiomorphic adenoma, a common type of benign epithelial tumor of the salivary gland. Multiple splice variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jun 2018]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).