

## Product datasheet for **TL311750V**

### galectin 9 (LGALS9) Human shRNA Lentiviral Particle (Locus ID 3965)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	galectin 9 (LGALS9) Human shRNA Lentiviral Particle (Locus ID 3965)
Locus ID:	3965
Synonyms:	HUAT; LGALS9A
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	LGALS9 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001330163</a> , <a href="#">NM_002308</a> , <a href="#">NM_009587</a> , <a href="#">NR_024043</a> , <a href="#">NM_002308.2</a> , <a href="#">NM_002308.3</a> , <a href="#">NM_009587.1</a> , <a href="#">NM_009587.2</a> , <a href="#">BC110340</a> , <a href="#">BC110340.1</a> , <a href="#">BC034392</a> , <a href="#">BC105942</a> , <a href="#">BC105944</a> , <a href="#">NM_009587.3</a> , <a href="#">NM_002308.4</a>
UniProt ID:	<a href="#">O00182</a>
Summary:	The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. The protein encoded by this gene is an S-type lectin. It is overexpressed in Hodgkin's disease tissue and might participate in the interaction between the H&RS cells with their surrounding cells and might thus play a role in the pathogenesis of this disease and/or its associated immunodeficiency. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).