

Product datasheet for **TL311760**

Leptin Receptor (LEPR) Human shRNA Plasmid Kit (Locus ID 3953)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Leptin Receptor (LEPR) Human shRNA Plasmid Kit (Locus ID 3953)
Locus ID:	3953
Synonyms:	CD295; LEP-R; LEPRD; OB-R; OBR
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	LEPR - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3953). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001003679 , NM_001003680 , NM_001198687 , NM_001198688 , NM_001198689 , NM_002303 , NM_001003680.1 , NM_001003680.2 , NM_001003680.3 , NM_001003679.1 , NM_001003679.2 , NM_001003679.3 , NM_002303.1 , NM_002303.2 , NM_002303.3 , NM_002303.4 , NM_002303.5 , NM_001198688.1 , NM_001198689.1 , NM_001198687.1 , BC031086 , BC035120 , BC131779 , BC160009 , NM_001198687.2 , NM_002303.6
UniProt ID:	P48357
Summary:	The protein encoded by this gene belongs to the gp130 family of cytokine receptors that are known to stimulate gene transcription via activation of cytosolic STAT proteins. This protein is a receptor for leptin (an adipocyte-specific hormone that regulates body weight), and is involved in the regulation of fat metabolism, as well as in a novel hematopoietic pathway that is required for normal lymphopoiesis. Mutations in this gene have been associated with obesity and pituitary dysfunction. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. It is noteworthy that this gene and LEPROT gene (GeneID:54741) share the same promoter and the first 2 exons, however, encode distinct proteins (PMID:9207021).[provided by RefSeq, Nov 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).