

## Product datasheet for **TL311917V**

### Kinesin Heavy Chain 2 (KIF2A) Human shRNA Lentiviral Particle (Locus ID 3796)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Kinesin Heavy Chain 2 (KIF2A) Human shRNA Lentiviral Particle (Locus ID 3796)
Locus ID:	3796
Synonyms:	CDCBM3; HK2; KIF2
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	KIF2A - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001098511</a> , <a href="#">NM_001243952</a> , <a href="#">NM_001243953</a> , <a href="#">NM_004520</a> , <a href="#">NM_004520.2</a> , <a href="#">NM_004520.3</a> , <a href="#">NM_004520.4</a> , <a href="#">NM_001098511.1</a> , <a href="#">NM_001098511.2</a> , <a href="#">NM_001243952.1</a> , <a href="#">NM_001243953.1</a> , <a href="#">BC031828</a> , <a href="#">BC047544</a> , <a href="#">NM_001098511.3</a>
UniProt ID:	<a href="#">O00139</a>
Summary:	The protein encoded by this gene is a plus end-directed motor required for normal mitotic progression. The encoded protein is required for normal spindle activity during mitosis and is necessary for normal brain development. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).