

Product datasheet for **TL312021**

KCNE1 Human shRNA Plasmid Kit (Locus ID 3753)

Product data:

Product Type:	shRNA Plasmids
Product Name:	KCNE1 Human shRNA Plasmid Kit (Locus ID 3753)
Locus ID:	3753
Synonyms:	ISK; JLNS; JLNS2; LQT2/5; LQT5; MinK
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	KCNE1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3753). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC036452 , NM_000219 , NM_001127668 , NM_001127669 , NM_001127670 , NM_001270402 , NM_001270403 , NM_001270404 , NM_001270405 , NM_000219.1 , NM_000219.2 , NM_000219.3 , NM_000219.4 , NM_001127668.1 , NM_001127668.2 , NM_001127668.3 , NM_001127669.1 , NM_001127669.2 , NM_001127669.3 , NM_001127670.1 , NM_001127670.2 , NM_001127670.3 , NM_001270402.1 , NM_001270402.2 , NM_001270403.1 , NM_001270403.2 , NM_001270404.1 , NM_001270404.2 , NM_001270405.1 , NM_001270405.2 , BC036452.1 , BC046224 , BC069055 , NM_000219.6 , NM_001127669.4
UniProt ID:	P15382
Summary:	The product of this gene belongs to the potassium channel KCNE family. Potassium ion channels are essential to many cellular functions and show a high degree of diversity, varying in their electrophysiologic and pharmacologic properties. This gene encodes a transmembrane protein known to associate with the product of the KVLQT1 gene to form the delayed rectifier potassium channel. Mutation in this gene are associated with both Jervell and Lange-Nielsen and Romano-Ward forms of long-QT syndrome. Alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008]



[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>