

## Product datasheet for **TL312080**

### Integrin beta 4 (ITGB4) Human shRNA Plasmid Kit (Locus ID 3691)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Integrin beta 4 (ITGB4) Human shRNA Plasmid Kit (Locus ID 3691)
Locus ID:	3691
Synonyms:	CD104; GP150
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ITGB4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3691). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_000213</a> , <a href="#">NM_001005619</a> , <a href="#">NM_001005731</a> , <a href="#">NM_001321123</a> , <a href="#">NM_001005731.1</a> , <a href="#">NM_001005731.2</a> , <a href="#">NM_000213.1</a> , <a href="#">NM_000213.2</a> , <a href="#">NM_000213.3</a> , <a href="#">NM_000213.4</a> , <a href="#">NM_001005619.1</a> , <a href="#">BC109238</a> , <a href="#">BC118916</a> , <a href="#">BC126411</a> , <a href="#">BC143738</a> , <a href="#">BC143742</a> , <a href="#">NM_001005731.3</a> , <a href="#">NM_000213.5</a>
UniProt ID:	<a href="#">P16144</a>
Summary:	Integrins are heterodimers comprised of alpha and beta subunits, that are noncovalently associated transmembrane glycoprotein receptors. Different combinations of alpha and beta polypeptides form complexes that vary in their ligand-binding specificities. Integrins mediate cell-matrix or cell-cell adhesion, and transduced signals that regulate gene expression and cell growth. This gene encodes the integrin beta 4 subunit, a receptor for the laminins. This subunit tends to associate with alpha 6 subunit and is likely to play a pivotal role in the biology of invasive carcinoma. Mutations in this gene are associated with epidermolysis bullosa with pyloric atresia. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).